


VITA In-Ceram® ALUMINA – Slip technique

All-ceramic system for fabricating highly stable crown and anterior three-unit bridge frameworks with one pontic



VITA

Product description

VITA In-Ceram ALUMINA ($\text{Al}_2\text{O}_3/\text{ZrO}_2$) is indicated for the fabrication of highly stable oxide-ceramic frameworks for both anterior and posterior single crowns and three-unit anterior bridges with one pontic.

After glass infiltration firing, during which the sintered frameworks retain their final strength and individual coloring, they are veneered with VITA VM 7 fine-structure ceramic.

VITA In-Ceram® ALUMINA kit



VITA In-Ceram ALUMINA kit Prod. no. HSORALV2

Quantity	Content	Material
1	25 g	ALUMINA GLASS POWDER AL2
1	25 g	ALUMINA GLASS POWDER AL4
1	400 g	ALUMINA POWDER
20	20 g	Special plaster for crowns / bridges
20	5 ml	ALUMINA / ZIRCONIA mixing liquid
1	30 ml	Interspace varnish
1	30 ml	Interspace varnish thinner
1	6 ml	Testing liquid
1	5 ml	ALUMINA additive
3	pcs	Duplicating rings, Ø 28 mm
1	pc	Vacuum mixing beaker, incl. rubber cap and glass tube
1	pc	Mixing beaker, acrylic
1	pc	Spare beaker with lid, acrylic
1	pc	Pipette, gradations 1/10 ml
1	pc	Glass spatula
1	pc	Brush no. IC 4
2	pcs	In-Ceram firing trays for bridge technique
1	pc	In-Ceram GLASS POWDER shade indicator
		Instructions for use

Printed materials

Title	Description	Order no.
VITA In-Ceram ALUMINA	Instructions for use	820

VITA In-Ceram® ALUMINA products

Description	Size	Content	Prod. no.
ALUMINA GLASS POWDER AL1	Bottle	25 g	HGAL125
ALUMINA GLASS POWDER AL2	Bottle	25 g	HGAL225
ALUMINA GLASS POWDER AL3	Bottle	25 g	HGAL325
ALUMINA GLASS POWDER AL4	Bottle	25 g	HGAL425
ALUMINA POWDER	Tin	200 g	HP200
ALUMINA POWDER	Tin	400 g	HP400
Special plaster for crowns / bridges	20 sachets each	20g	HGN20
Interspace varnish	Bottle	30 ml	HD30
Interspace varnish thinner	Bottle	30 ml	HDV30
ALUMINA / ZIRCONIA mixing liquid	20 ampoules each	5 ml	HAFN5
ALUMINA additive	Bottle	5 ml	HA5
Testing liquid	Bottle	6 ml	HP6
Duplicating ring	Ø 28 mm	1 pc	B019
Glass spatula	-	1 pc	B056
Pipette	Gradations 1/10 ml	1 pc	B052
Brush no. IC 4	-	1 pc	B092
Mixing beaker with lid, acrylic	100 ml	1 pc	B003
Spare beaker with lid, acrylic	-	1 pc	B099
Vacuum mixing beaker, incl. rubber cap and glass tube	-	1 pc	B130
Firing tray for bridge technique	pk	2 pcs	B007
In-Ceram GLASS POWDER shade indicator	-	1 pc	B2711C

Features/Benefits

- A highly translucent all-ceramic system which has clinically proven itself millions of times over since 1989.
- Due to their high strength of around 500 MPa, restorations made from VITA In-Ceram ALUMINA offer the possibility of non-adhesive cementation.